

POCKET WOOD LIGHT



|               |                |
|---------------|----------------|
| Product Code  | 31189          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|         |               |
|---------|---------------|
| EAN/UPC | 8023279311891 |
| PARAF   | 903625731     |

Description

Pocket ultraviolet light, to detect skin disorders.  
Provided with replaceable fluorescent lamp.  
Power 6 V.  
Works with 4 x AA batteries with 4-6 hours autonomy.  
Size: 16 x 5 x 2 cm.