


VISTA WOOD LED LAMP - handle



Product Code	31195
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CND	K02010202

Description

VISTA 4 LED WOOD LAMP 3W - handle

Lamps to diagnose the condition of skin that appears in different colours under ultraviolet rays.
Polished glass magnifying lens 10x6 cm - 3 dioptries.

LED technology emits two types of light: wood's light at 365 nm wavelength (UV) and daylight (CRI >80%, 4,500°K) that can be switched by button.

Voltage: 100-240 V - 50/60 Hz
Manual in GB, FR, IT, ES, PT.

Handle size: Ø 3x12 cm
Upper side with protective sliding cover
Lower side with 4 LEDs (2 UV + 2 white)

Made in Spain.