



|               |                |
|---------------|----------------|
| Product Code  | 31197          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |               |
|------------|---------------|
| UK-REP     | No            |
| CH-REP     | No            |
| RDM (NSIS) | 1887267       |
| CND        | Z120614       |
| EAN/UPC    | 8436562860028 |

## Description

VISTA WOOD LAMP 9W - trolley

Lamps to diagnose the condition of skin that appears in different colours under ultraviolet rays. Polished glass magnifying lens 10x6 cm - 3 diopres.

Fluorescent lamps have 9W special compact fluorescent tube at approximately 360 nm wavelength.

Voltage: 220-240 V - 50/60 Hz  
Manual in GB, FR, IT, ES, PT.

Arm length 45+45 cm  
Trolley height 95 cm

Made in Spain.