

NYSTAGMUS SPECTACLES FRENZEL



Product Code	31204
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2398396
CND	Z12120112
EAN/UPC	8023279312041
PARAF	903145757

Description

NYSTAGMUS SPECTACLES - DR. FRENZEL - black

Easy to use in the search of spontaneous or provoked nystagmus.
Fixed, biconvex lens + 18.0 diopters.
Nystagmus spectacles should be connected only with < 4.2 V transformer as per safety norms 93/42/EEC. See code 31203 or 31214.

Suggested for long examinations
Supplied with head strip
Multilingual manual: GB, FR, IT, ES, PT, GR

Made in Germany.