

LED BATTERY OPERATED NYSTAGMUS SPECTACLES

- white



Product Code	31207
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2401556
CND	Z12120112
EAN/UPC	8023279312072

Description

LED BATTERY NYSTAGMUS SPECTACLES - white

With LED and inner function that makes easier watching the rashes within a particular time period (30, 60, 90 seconds).

This model allows wireless working and is provided with power through the adjustable battery placed above the spectacle. Requires three 1.5 Volt batteries, size C LR14/AM-2 (baby cells). Illumination can be adjusted through a turn of the switch.

Suggested for short examinations
Adjustable illumination through a slight turn around of the switch
Multilingual manual: GB, FR, IT, ES, PT, GR