

VIDEO-NYSTAGMOSCOPE - USB - mask with Ir camera



Product Code	31213
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2190611
CND	Z12140401
EAN/UPC	8023279312133

Description

VIDEO-NYSTAGMOSCOPE - USB - mask with Ir camera and USB cable

Facilitate to view and register the Nystagmus through standard vestibular tests (both position and instrumental proof). This system allows to view on the monitor or PC/Notebook images of patient eyes during tests and to make Visual Suppression Test (VST) through a green LED light inside the mask to allow fixation. Innovative technology allows to get quality images. Power 12V DC.

752x582 pixelsresolution 520 TV lines

Made in Italy.

752x582 pixels resolution 520 TV lines.

Technical Specifications

752 x 582 pixel resolution; 520 TV lines