


ALUMINIUM TUNING FORK 128 Hz



Product Code	31225
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2533755
CND	V030208
EAN/UPC	8023279312256
PARAF	964181844
GTIN	08023279312256
GMDN	32526

Description

ALUMINIUM TUNING FORK 128 Hz

Made of finest aluminium alloy and precisely shaped to meet highest requisites (standards) in medical field.

20.4 cm long