


ALUMINIUM TUNING FORK 128 Hz - with foot



Product Code	31226
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2533755
CND	V030208
EAN/UPC	8023279312263
GTIN	08023279312263
GMDN	32526

Description

ALUMINIUM TUNING FORK 128 Hz - with foot

Made of finest aluminium alloy and precisely shaped to meet highest requisites (standards) in medical field.

- Resting pad for forehead
- 22 cm long