

DIAPASON HARTMANN 256 Hz



Product Code	31231
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2733113
CND	V030208
EAN/UPC	8023279312317
UMDNS	14-255
PARAF	909921951
GTIN	08023279312317
GMDN	32526

Description

ENT TUNING FORKS

- Tuning fork 256 Hz - with weight

All made in Germany, with age-long craftsmanship experience, in special alloy that ensures maximum sound precision.

Made in Germany.