

SET DIAPASON HARTMANN - 5 pcs.



Product Code	31240
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2733113
CND	V030208
EAN/UPC	8023279312409
UMDNS	14-255
PARAF	903120196
GTIN	08023279312409
GMDN	32526

Description

ENT TUNING FORKS

- Set of 5 tuning forks in a stylish lined case (contains 128 Hz*, 256 Hz*, 512 Hz, 1024 Hz, 2048 Hz).
Tuning forks with asterisk * are supplied with weight

All made in Germany, with age-long craftsmanship experience, in special alloy that ensures maximum sound precision.

Made in Germany.