

UNISEX COTTON SOCKS - M - strong compression - black



| | |
|---------------|----------------|
| Product Code | 41241 |
| Unit of sale | 1 pair |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| UK-REP | No |
| CH-REP | No |
| RDM (NSIS) | 2636014 |
| CND | M030405 |
| EAN/UPC | 8033201402658 |
| PARAF | 925565095 |

Description

UNISEX GRADUAL COMPRESSION COTTON SOCKS - pair

- Strong compression socks
- Colour: Black
- Size: Medium
- Shoe size: [European 39-41] [US Men: 6½-8] [US Women: 8-10]

High quality knee socks made of soft cotton combined with elastic fibres to provide proper gradual compression.

Special attention was applied to the development of this cotton product line to ensure maximum fibre quality, along with a soft elastic cuff that offers excellent levels of expandibility: all to ensure maximum comfort and wellness throughout the day.

This product line is available in two compression levels (medium and strong) and in two colours (black and blue). On request, for quantities of 100 or more they are also available in light and extra strong compression and beige and grey colours.

Medium compression level: 14-16 mm/Hg at the ankle

To help prevent discomfort and heaviness in the legs. Composition: 67% cotton, 26% polyamide, 7% elastane.

High compression level: 15-21 mm/Hg at the ankle

To help prevent the onset of venous disease. Composition: 63% cotton, 27% polyamide, 10% elastane.

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, GR, HU

Made in Italy.