

WORTH LED LIGHT TEST with red/green glasses



Product Code	31289
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

EAN/UPC	8906200830602
GMDN	12293

Description

WORTH LED LIGHT TEST - with red/green glasses

The Worth 4 dot test, is a clinical test for suppression of either the right or left eye. Suppression occurs during binocular vision, when the brain does not process the information received from either of the eyes. This is a common adaptation to strabismus, amblyopia, and aniseikonia.
Battery not included, glasses included.

Multilingual instructions: GB, FR, IT, ES, PT, DE, GR, Arabic.