



Product Code	31300
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2753418
CND	Z12149006
EAN/UPC	8023279313000
PARAF	903145819
GTIN	08023279313000
GMDN	12849

Description

Dual function instrument features convertible diagnostic pen light and mini-otoscope with 3 different sizes 2.5 - 3.5 - 4.5 x 38 mm of funnels or specula.

White LED light. Packed in a gift box.

Pupil gauge and dual slide rule printed on the penlight.

Supplied with batteries (2 x AAA).