



| | |
|---------------|-------------------|
| Product Code | 21405 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|---------|---------------|
| EAN/UPC | 8023279214055 |
| PARAF | 977809627 |

Description

DOCTOR'S WHITE COATS 60% COTTON - 40% POLYESTER

- WOMAN - With 5 buttons
- Size: (Europe: 54-56) (USA: XL) (France: 50-52)
- 60% cotton / 40% polyester

Professional and stylish line with all sizes for men and women, made of 60% high quality cotton and 40% polyester.

Doctor's white coats have 2 side pockets and 1 breast pocket.

Man and woman model have buttons, unisex model has hidden press studs and elasticated cuffs.

Weight 200 g/m²

Technical Specifications

Composition: 60% cotton - 40% polyester

Weight: 200 g/m²

Yarn count: warp: Ne 20 ± 2 / weft: Ne 18 ± 2

Construction: warp: 40 ± 3 / weft: 19 ± 3

Breaking strength - warp: 100 kg ±3%

Breaking strength - weft: 39 kg ±3%

Fastness to washing: 4

Treatments and finishing: crease resist

Dimensional stability - washing at 90°C: 3%