

# MULTIFUNCTIONAL OTOSCOPE 4 in 1 - manual in 24 languages



Product Code	31410
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2804444
CND	Z12021009
EAN/UPC	6976313750809
PARAF	950850659

## Description

4 IN 1 MULTIFUNCTIONAL OTOSCOPE - ear, nose, throat, mouth

4 in 1 handheld otoscope equipped with lighting and a camera for inspection of ear, nose, throat and mouth. Users can capture high-resolution images and videos to better examine and document their findings for easy sharing. Its wireless connectivity enables real-time viewing and sharing of examinations through Deepview Otoscope App. The built-in LED light provides sufficient illumination to ensure clear view. The device is suitable for children over 1 year old and adults.

- 0.3 MP HD camera
- 4.0 mm ultra slim lens
- 3-Axis Gyroscope
- 360° view
- Waterproof lens

Supplied with 2 silicone nose specula, 2 silicone ear specula, 1 tongue depressor, charging cable and operating manual in GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic.

## Technical Specifications

Operating frequency: 2.4 GHz  
Image transmission rate: 30 fps  
LEDs: cold 6 LED lights  
Camera diameter: 4.0 mm  
Optimum focal length: 18 mm  $\pm$ 5 mm  
Pixels: 0.3 M  
Power input: DC 5V 1A  
Power supply: 800 mAh lithium battery  
Battery life: 4 hours  
Size - weight: 23x57x155 mm - 60 g  
Protection: internal power supply, type BF. IP22