

WHITE COAT WITH STUD - cotton/polyester - unisex size XXXL



Product Code	21427
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

EAN/UPC	8023279214277
PARAF	977809755

Description

DOCTOR'S WHITE COATS 60% COTTON - 40% POLYESTER

- UNISEX - With press stud buttons
- Size: (Europe: 62-64) (USA: XXXL) (France: 58-60)
- 60% cotton / 40% polyester

Professional and stylish line with all sizes for men and women, made of 60% high quality cotton and 40% polyester.
Doctor's white coats have 2 side pockets and 1 breast pocket.
Man and woman model have buttons, unisex model has hidden press studs and elasticated cuffs.
Weight 200 g/m²

Technical Specifications

Composition: 60% cotton - 40% polyester
Weight: 200 g/m²
Yarn count: warp: Ne 20 ± 2 / weft: Ne 18 ± 2
Construction: warp: 40 ± 3 / weft: 19 ± 3
Breaking strength - warp: 100 kg ±3%
Breaking strength - weft: 39 kg ±3%
Fastness to washing: 4
Treatments and finishing: crease resist
Dimensional stability - washing at 90°C: 3%