

SET 3 VETERINARY SPECULUM



Product Code	31447
Unit of sale	set of 3 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279314472

Description

AUTOCLAVABLE SPECULUM FOR OTOSCOPES - LATEX FREE

Kit 3 autoclavable veterinary speculum Ø 4, 5, 7 mm

To be used with GIMA code: 31442

- made from resilient thermoplastic material
- smooth matte finish surface ensure easy cleaning
- metallic bayonet lock for positive fit to otoscope
- autoclavable at 134°C