

SET 3 ADULT SPECULUM



| | |
|---------------|----------------|
| Product Code | 31448 |
| Unit of sale | set of 3 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2831169 |
| CND | Q030499 |
| EAN/UPC | 8023279314489 |
| GTIN | 08023279314489 |
| GMDN | 33395 |

Description

AUTOCLAVABLE SPECULUM FOR OTOSCOPES - LATEX FREE

Kit 3 autoclavable speculum Ø 2.5, 3.5, 4.5 mm

To be used with GIMA code: 31435-37, 31440-42, 31471, 31473, 31475, 31477, 31480
Also compatible with: Heine (Non compatibile con Mini 3000), Kawe, Riester

- made from resilient thermoplastic material
- smooth matte finish surface ensure easy cleaning
- metallic bayonet lock for positive fit to otoscope
- autoclavable at 134°C