


LIGHTED MIXED CURETTES with lens and LED illuminator



Product Code	31450
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	215023
CND	L149004
EAN/UPC	8023279314502
GMDN	65446

Description

KIT MIXED LIGHTED CURETTES

Kit is composed by 50 lighted currettes (8 FlexLoop and 7 of other 6 tip styles) + LED illuminator + magnification lens