



Product Code	31460
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2681599
CND	Q030499
EAN/UPC	0014646012228
PARAF	923363612
GMDN	65445

Description

MICROLOOP cures: For compact cerumen, more rigid construction, Tip can be angled for added strength. Rounded tip for comfort 3 mm wide tip.

Our four styles of cures give you the most effective, economical and safe way to remove cerumen - whatever the size of ear or consistency of cerumen.

- Patient comfort: metal cures are cold, hard and sharp. Blonlx ear cures with rounded plastic tips reduce trauma and risk of laceration. Latex free.
- Effectiveness: ear cures can duplicate the effectiveness of metal cures in almost any situation and their polypropylene composition makes them virtually unbreakable.
- Disposability: ear cures don't require expensive and time-consuming sterilization and eliminate the risk of cross-contamination.