

# XENON-HALOGEN OPERATIVE DIAGNOSTIC SET - 3.5

V



Product Code	31482
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2831176
CND	Z121402
EAN/UPC	8023279314823
GTIN	08023279314823
GMDN	45084

## Description

### XENON-HALOGEN OPERATIVE DIAGNOSTIC SET - 3.5 V

Supplied in plastic rigid case containing: one rechargeable handle, one otoscope head operative, one ophthalmoscope head, 3 autoclavable speculum Ø 2.5, 3.5, 4.5 mm, two laryngeal mirrors (size 3 and 4), adjustable nasal speculum, bent throat lamp, tongue spatula. Should be used with 31541 + 31542.