

INFLATION BULB FOR OTOSCOPES - grey



Product Code	31483
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2324878
CND	Q0399
EAN/UPC	8023279314830
GTIN	08023279314830

Description

Pneumatic PVC bulb for Oscopes

To be used with GIMA code: 31435/36/37/40/42/71/73/75/77-31520/21/22, 31535

Also compatible with: Heine, Kawe, Riester

Weight 60g

Total length 54 cm, tube 46 cm