

EAR SPECULUM diam. 4 mm for WelchAllyn - disposable - black



Product Code	31509
Unit of sale	box of 250 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2078999
CND	Q030499
EAN/UPC	3701126705318
GTIN	03701126705318

Description

Disposable ear-speculum Ø 4 mm - black

To be used with:
Welch Allyn Pocket-Pro and Pocket scope Series 200, 250

Multilingual box: GB, FR, IT, ES, PT, DE, GR, DK, Arabic.