

# SIGMA F.O. OPHTHALMOSCOPE 2.5V - Xenon-halogen - case



Product Code	31527
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2831183
CND	Z12120114
EAN/UPC	8023279315271
UMDNS	46788
GTIN	08023279315271
GMDN	45084

## Description

- SIGMA F.O. OPHTHALMOSCOPE 2.5V - standard handle
- Black
- Xenon-Halogen
- Packing: Case

High quality pocket diagnostic instruments designed for professional use in ENT and ophthalmology. Light weight and manufactured from impact resistant re-enforced plastic.

Manufactured with high quality materials and top-notch workmanship ensures excellent reliability in every examination.

- distal fibre optic illumination with super bright Xenon/Halogen light
- swivel type window with 3X magnification lens for sharp images and minimal reflection
- insufflation port for pneumatic testing of tympanic mobility