



Product Code	31529
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2831183
CND	Z12120114
EAN/UPC	8023279315295
GTIN	08023279315295
GMDN	45084

Description

- SIGMA F.O. OPHTHALMOSCOPE 2.5V - standard handle
- Black
- LED
- Packing: Case

High quality pocket diagnostic instruments designed for professional use in ENT and ophthalmology. Light weight and manufactured from impact resistant re-enforced plastic.

Manufactured with high quality materials and top-notch workmanship ensures excellent reliability in every examination.

Battery life with Led bulb 6 time longer than with traditional bulb

- distal fibre optic illumination with super bright LED light
- swivel type window with 3X magnification lens for sharp images and minimal reflection
- insufflation port for pneumatic testing of tympanic mobility

Technical Specifications

LED FEATURES

- colour temperature: 5,000°K
- CRI > 80
- real colours definition
- 85,000 LUX
- bulb estimated life 50,000 hours