



Product Code	31539
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2734804
CND	L1699
EAN/UPC	8023279315394
UMDNS	12849 + 46788
GTIN	08023279315394
GMDN	45084

Description

VISIO 2000 F.O. XENON DIAGNOSTIC SET - 3.5 V

F.O. otoscope head, F.O. ophthalmo head, battery handle with bayonet lock, nasal illuminator with F.O. light, curved illuminator with bulb, plain mirror 20.0 mm, plastic tongue depressor, 1 set of 5 reusable eartips (2.5, 3.0, 4.0, 5.0 & 10.0 mm) tube of 14 disposable eartips (2.5 & 4.0 mm), spare oto and ophthalmo bulbs, in a rigid case.

Professional diagnostic sets with fiber optic illumination, Xenon light 3.5 V, rheostat to adjust the intensity of light. Made of resistant re-inforced plastic and durable easy locking system. Supplied in deluxe case, should be used with charging station 31542 and Li-ion rechargeable battery 31541.

Technical Specifications

Xenon otoscope

- Swivel type window with 3X magnification for sharp images and minimal reflection
- Distal F.O. illumination with super bright LED light
- Durable easy locking system
- Rheostat

Xenon ophthalmoscope

- 6 apertures to select from:
 - Micro spot aperture
 - Small aperture
 - Large aperture
 - Fixation aperture
 - Slit
 - Red-free filter
- Enhanced Xenon super bright light
- Durable easy locking system
- Rheostat