

# E-SCOPE OTOSCOPE - halogen 2.7V - black in pouch



|               |                |
|---------------|----------------|
| Product Code  | 31548          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|            |               |
|------------|---------------|
| RDM (NSIS) | 1830774       |
| CND        | Z12149006     |
| EAN/UPC    | 4045396176898 |

## Description

### E-SCOPE® OTOSCOPE - HALOGEN 2.7V - BLACK IN POUCH

Halogen illumination otoscope.  
Sealing precision optics pivoting to both sides with triple magnification for superior view.  
Suitable for performance of pneumatic tests (supplied without connector and insufflator).