

E-SCOPE DIAGNOSTIC KIT - LED 3.7V - white in rigid case



Product Code	31558
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

EAN/UPC	4045396177031
---------	---------------

Description

E-SCOPE® DIAGNOSTIC KIT - WHITE
Oto-ophthalmoscope sets, LED 3.7V, in rigid case.

- **Otoscope**
Sealing precision optics pivoting to both sides with triple magnification for superior view.
Suitable for performance of pneumatic tests (supplied without connector and insufflator).
- **Ophthalmoscope**
 - innovative 3.7V LED illumination for higher-contrast and more efficient diagnosis
 - dioptre disc with 18 corrective lenses ($\pm 1, 2, 3, 4, 6, 8, 10, 15, 20$)
 - easy-to-use aperture wheel with six apertures (fixation star, large circle, small circle, red free filter, blue filter, semi circle)
 - high performance optics with aspherical condenser
 - spectacle protection
 - parallel optical path

INNOVATIVE LED 3.7 V ILLUMINATION :

- real colour 5,500K white light
- long life LED bulbs: up to 20,000 hours
- low battery consumption: up to 50 hours operation through IPC (Integrated Power Converter). It converts the voltage of AA batteries into the optimum voltage required for LED