

SIGMA F.O. LED OTO-OPHTHALMOSCOPE SET with 2 rechargeable handles - case



| | |
|---------------|----------------|
| Product Code | 31584 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2734775 |
| CND | L1699 |
| EAN/UPC | 8023279315844 |
| GTIN | 08023279315844 |
| GMDN | 45084 |

Description

Oto-ophthalmoscope set with 2 rechargeable handles and batteries

Type: LED

Packing: Case

This set is supplied with oto and ophthalmo head, 2 handles, 14 disposable specula (7 with Ø 2.5 mm, 7 with Ø 4 mm), 2 rechargeable batteries.