



Product Code	21605
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

EAN/UPC	8023279216059
---------	---------------

## Description

JACKET WITH STUDS 60% COTTON - 40% POLYESTER

- Size: XL
- Colour: white
- UNISEX
- Five hidden press studs
- Short sleeves
- 2 front and 1 breast pockets

Professional and stylish line with all sizes for men and women, made of 60% high quality cotton and 40% polyester.  
Stock in 5 colours, other 22 colours available in 60 days.  
White jacket 200 g/m<sup>2</sup>, coloured jacket 185 g/m<sup>2</sup>.

## Technical Specifications

Composition: 60% cotton - 40% polyester  
Weight: 200 g/m<sup>2</sup> (white), 185 g/m<sup>2</sup> (coloured)  
Yarn count: warp: Ne 20 ± 2 weft: Ne 18 ± 2  
Construction: warp: 40 ± 3 weft: 19 ± 3  
Breaking strength - warp: 100 kg ±3%  
Breaking strength - weft: 39 kg ±3% (white) 49 kg ±3% (coloured)  
Fastness to washing: 4  
Treatments and finishing: crease resist  
Dimensional stability - washing at 90°C: 3%