

HEINE MINI 3000 LED DIRECT OTOSCOPE with
reusable tips and case - blue - D-856.20.021



Product Code	31665
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

EAN/UPC	4053755202645
---------	---------------

Description

HEINE MINI 3000 LED 2.5V OTOSCOPE

Mini 3000 LEDHQ direct illumination otoscope

- Blue
- Code HEINE: D-856.20.021
- Reusable Tips 4 pcs x 4 mm
- AllSpec disposable Tips 10 pcs x 4 mm
- Battery handle (Batteries included)
- Hard case

Compact otoscope with direct LED illumination. Natural colour rendering, no need to ever change the LED. Exclusive design.

Detailed view

- it eliminates reflections due to the head's black inner coating
- high-resolution images
- colour temperature 4,000°K
- 3x magnification
- particularly bright due to LEDHQ

Intuitive operability and practical attachment clip

- it switches off automatically as soon as it is clipped to the pocket
- easy instrumentation. Viewing window fixed , can be swivelled to the side
- operating time up to 10 hours