

# HEINE MINI 3000 F.O. LED OTOSCOPE with reusable tips and case - blue - D-885.20.021



Product Code	31679
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

EAN/UPC 4053755161188

## Description

### HEINE MINI 3000 LED FIBER OPTICS (F.O.) OTOSCOPES

#### Mini 3000 LEDHQ F.O. otoscope

- Blue
- Code HEINE: D-885.20.021
- Reusable Tips 4 pcs x 4 mm
- AllSpec disposable Tips 10 pcs x 4 mm
- Battery handle (Batteries included)
- Hard case

Compact otoscope with especially bright illumination and authentic colour rendering due to LEDHQ (no need to ever change the LED). Exclusive design and sturdy housing.

#### Detailed view

- high-resolution images, minimal reflection
- 3x magnification
- reflex-free illumination of tympanum and ear canal due to fiber optic technology
- particularly bright due to LEDHQ
- operating time up to 10 hours

#### Intuitive operability and practical attachment clip

- it switches off automatically as soon as it is clipped to the pocket
- easy instrumentation. Viewing window fixed, can be swivelled to the side
- pneumatic eardrum testing by means of an integrated insufflation port