

HEINE MINI 3000 F.O.OTOSCOPE - rechargeable handle
- black



Product Code	31682
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	97394
CND	Z12149006
EAN/UPC	4053755178704
UMDNS	12-849
GMDN	12849

Description

MINI 3000 F.O. XENON OTOSCOPES - RECHARGEABLE HANDLE

- F.O. Xenon
- RECHARGEABLE HANDLE (Need Heine Mini 3000 NT charger code 31692)
- Colour: Black
- HEINE code: D-001.70.301
- Supplied in a box with 5 disposable speculum Ø 4 mm and 5 Ø 2.5 mm

Powerful pocket otoscope in a unique compact design.

- Xenon halogen bulbs for bright illumination and accurate colour rendering
- 20,000 on/off switch cycles guaranteed
- distal Fiber Optic (F.O.) Illumination: Reflex-free illumination of tympanum and ear canal
- viewing window with 3X magnification and optimized casing surface for razor-sharp images and minimal reflection. Swivelling viewing window built into instrument
- pneumatic testing of the tympanum mobility by using the insufflation port and an optional connector
- attachment clip with integrated on/off switch
- replaceable batteries: size AA or optional rechargeable battery (31693) with mini NT table charger (31692).
- 2-piece handle and head. Easy to maintain, flexible. Compatible with other handles.