



Product Code	31688
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	Yes
EAN/UPC	4053755158737

Description

HEINE MINI 3000 LED 2.5V OPHTHALMOSCOPES

- Mini 3000 LEDHQ ophthalmoscope
- Black
 - Code HEINE: D-008.71.120
 - Battery handle (Batteries included)

Modern pocket ophthalmoscope with LEDHQ. Concentrated, bright light for optimum illumination, a brilliant image and authentic colour rendering. Fixation star aperture with grey filter to reduce light intensity. Excellent optics and outstanding colour rendering with the smallest dimensions.

- Reliable diagnosis
- simple examination process possible without pupil dilation
 - clear, high-resolution image with minimal reflections
 - 18 lenses ranging from - 20 to + 20 D prevent higher refractive errors
 - fixation star including grey filter for patients who are sensitive to light
 - 5 different apertures

- Intuitive operability
- it switches off automatically as soon as it is clipped to the pocket
 - practical one-finger operation
 - handle and head are in 2 parts for easy servicing and expandability