



Product Code	41700
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2001716
CND	Z12080303
EAN/UPC	6970392218809
UMDNS	10485
PARAF	984594895

Description

MANUAL BREAST PUMP

Compact and lightweight breast pumps for efficient and discreet pumping. Easy to assemble, use and clean. Supplied with breast shield, nipple with cover, milk bottle with lid, multilingual box and manual (GB, FR, IT, ES, PT, DE).

Technical Specifications

Vacuum: - 250 mmHg
Weight: 159 g
Size: 205 x 186 x 90 mm
Operation: manual