

SINGLE ELECTRIC BREAST PUMP



Product Code	41701
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2001713
CND	Z12080303
EAN/UPC	6970392218908
UMDNS	10485
PARAF	984594883

Description

SINGLE ELECTRIC BREAST PUMP

Compact and lightweight breast pumps for efficient and discreet pumping. Easy to assemble, use and clean. Supplied with breast shield, nipple with cover, milk bottle with lid, multilingual box and manual (GB, FR, IT, ES, PT, DE).

All parts that touch breast milk are made without BPA, phthalate free.

2 phase design, easy for stimulation and expression.

10 vacuum levels

Soft and comfortable breast shield

Graduated bottle (ml/Oz)

Technical Specifications

Vacuum Stimulate: - 60-130 mmHg

Vacuum Expression: - 60-250 mmHg

Weight: 265 g

Size: 129 x 129 x 55 mm

Operation: 4 AA batteries not included or AC adapter (optional)