

HEINE BETA 200 LED F.O. OTOSCOPE - B-141.10.118



Product Code	31724
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

CH-REP	Yes
RDM (NSIS)	2780678
CND	Z12149006
EAN/UPC	4053755115402

Description

HEINE BETA 200 LEDHQ F.O. OTOSCOPE - 2.5V - B-141.10.118

Exceptional durability combined with optimum 3x magnification.
Particularly bright illumination thanks to LEDHQ.
Solid metal housing.

Supplied in hard case with BETA battery 2.5V handle with rheostat, a set of single use specula and a set of reusable specula.

Made in Germany.