

HEINE BETA 400 LED F.O. OTOSCOPE - B-143.24.420



Product Code	31728
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	2779104
CND	Z12149006
EAN/UPC	4053755182718
UMDNS	12-849
GMDN	12849

Description

HEINE BETA 400 LED F.O. OTOSCOPE - B-143.24.420

Exclusive 4.2X magnification.
LED in HQ with continuous brightness control.
All metal housing and scratch resistant glass lenses.
Colour temperature: 3,500 K. Colour rendering index (CRI) >97.
Exceptional depth of field, provides sharp spatial image.
LED otoscope with rechargeable handle supplied with 4 reusable speculum, 5 of both All Spec disposable speculum 2.5 and 4 mm Ø in a case with NT4 table charger (34374).

Made in Germany.