

**HEINE BETA 200 F.O. OPHTHALMOSCOPE - 3,5V Li-Ion
- C-144.27.388**



Product Code	31731
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	1725656
CND	Z12019080
EAN/UPC	4053755184422
UMDNS	12-817

Description

HEINE BETA 200 F.O. OPHTHALMOSCOPE - 3.5 V Li-ion - C-144.27.388

With USB rechargeable handle for mains socket and Li-ion batteries.

The unique aspherical optic system eliminates corneal and iris reflexes. The illuminated area of the retina remains fully visible even in case of small pupil (bright light). BETA 200 has 6 diaphragms with separately "green filter".
It is recommended for small or expanded pupils; suitable for slit or fluorescent ophthalmoscopy as well as for fixation test.

Made in Germany.
Range -35D to +40D. Supplied in hard case with replacement bulb.