

# HEINE BETA 200 F.O. OPHTHALMOSCOPE - 3,5V Li-Ion - C-144.28.388



|               |                |
|---------------|----------------|
| Product Code  | 31732          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|         |               |
|---------|---------------|
| CND     | Z12120114     |
| EAN/UPC | 4053755184514 |
| UMDNS   | 12-817        |
| GMDN    | 46786         |

## Description

HEINE BETA 200 F.O. Ophthalmoscope - 3.5 V Li-ion - C-144.28.388  
With USB rechargeable handle for mains socket and Li-ion batteries.

The unique aspherical optic system eliminates corneal and iris reflexes.

The illuminated area of the retina remains fully visible even in case of small pupil (bright light).

BETA 200 has 6 diaphragms with separately "green filter". It is recommended for small or expanded pupils; suitable for slit or fluorescent ophthalmoscopy as well as for fixation test.

Range -35D to +40D.

Supplied in hard case.

6 standard diaphragms: slit, fixation star with polar co-ordinates, cobalt blue filter, large spot, small spot, hemi spot

Made in Germany.