

# HEINE BETA 200 F.O. OPHTHALMOSCOPE SET - 3.5V Li-Ion C-144.27.376



Product Code	31736
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	104388
CND	Z12120114
EAN/UPC	4053755121915
UMDNS	12-817
GMDN	10284

## Description

Same as 31730 but with rechargeable handle and Li-ion batteries.

## Technical Specifications

- 6 standard diaphragms: slit, fixation star with polar co-ordinates, cobalt blue filter, large spot, small spot, hemi spot.
- 6 diaphragms with green filter.