

HEINE BETA 200S LED OPHTHALMOSCOPE - 3.5V C-261.28.388



Product Code	31737
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	1725637
CND	Z12019080
EAN/UPC	4053755184330
UMDNS	12-817

Description

HEINE BETA 200S LED OPHTHALMO - 3.5 V with USB rechargeable handle - C-261.28.388

- Superior aspherical optics and 74 single-diopter steps
- optical system: corneal and iris reflexes are eliminated. A full view is guaranteed even with small pupils.
 - XHL (Xenon Halogen Technology) or LED: bright, white light.
 - recessed, multi-coated viewing window: avoids stray light.
 - optics mounted on precision metal chassis.
 - ergonomic shape: instrument fits the orbita comfortably in any position.
 - 7 apertures with separate red-free filter:
 - lens range: perfect focus even with high refractive errors.
 - exact colour rendering
- Made in Germany.