

HEINE K 180 LED OPHTHALMOSCOPE - 2.5V



Product Code	31742
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

CH-REP	Yes
--------	-----

RDM (NSIS)	2779662
CND	Z12120114
EAN/UPC	4053755122356

Description

HEINE K180 LED OPHTHALMOSCOPE - 2.5 V - C-182.10.118

27 lenses from -35D to +40D
5 different apertures in 2 variants.

Made in Germany.