

HEINE BETA 200 F.O. LED RETINOSCOPE HEAD - 3.5 V  
- C-008.15.353



|               |                |
|---------------|----------------|
| Product Code  | 31748          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|         |               |
|---------|---------------|
| CND     | Z12120118     |
| EAN/UPC | 4053755191598 |

Description

HEINE code: C-008.15.353  
RETINOSCOPES HEADS: BETA 200 LED retinoscope head - 3.5 V  
Compatibility: Rechargeable handles