

RI-ACCU L Li-Ion BATTERY for C size - 10691



Product Code	31841
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2392866
CND	Z129004
EAN/UPC	4045396141506

Description

RI-ACCU L 3.5 V Li-ion BATTERY for "C" size - 10691

Li-ion batteries for great power
Higher charging capacity, longer lifetime, light-weight are the big advantages of Li-ions batteries. High performance for up to 1200 charging cycles. LEDs powered by lithium batteries can burn for up to 50 hours.

Made in Germany.