



Product Code	21853
Unit of sale	box of 100 pcs
Minimum order	1
Type	PPE
Class	III

UK-REP	Not required
EAN/UPC	5906615011811
PARAF	987476429

## Description

SANTEX® POWDERED - ambidextrous

Large size - non sterile  
Creamy colour

Non-sterile, powdered, examination and protective glove for single use.  
Made of natural rubber latex.  
Smooth external surface guarantees a tactile sensitivity during the work.  
They minimize hand sweating.  
Enhanced tactile sensitivity.  
AQL=1.0

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, GR, SE, NO, FI, DK, SK, CZ, BG, RO, HU, HR, RS, SL, KZ, LT, LV, RU, TR.

## Technical Specifications

External surface: Smooth  
Cuff: beaded  
MD (Medical Device): Class I  
PPE (Personal protective equipment): Category III