

SANTEX LATEX GLOVES - PRE POWDERED - extra large



Product Code	21854
Unit of sale	box of 100 pcs
Minimum order	1
Type	PPE
Class	III

UK-REP	Not required
EAN/UPC	5906615011828
PARAF	950454658

Description

SANTEX® POWDERED - ambidextrous

Extra large size - non sterile
Creamy colour

Non-sterile, powdered, examination and protective glove for single use.
Made of natural rubber latex.
Smooth external surface guarantees a tactile sensitivity during the work.
They minimize hand sweating.
Enhanced tactile sensitivity.
AQL=1.0

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, GR, SE, NO, FI, DK, SK, CZ, BG, RO, HU, HR, RS, SL, KZ, LT, LV, RU, TR.

Technical Specifications

External surface: Smooth
Cuff: beaded
MD (Medical Device): Class I
PPE (Personal protective equipment): Category III