

RIESTER ELITEVUE LED OTOSCOPE 3.5V - 2200-203



Product Code	31860
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC	4045396209008
PARAF	950833943

Description

RIESTER ELITEVUE LED OTOSCOPE 3.5 V - 2200-203

Both with C-handle for rechargeable Li-Ion battery
EliteVue offers twice the field of view and 66% greater magnification (at 15 mm distance) compared to standard otoscopes.
Improved view of auditory canal and tympanic membrane guarantees more efficient diagnosis.
Suitable for throat illumination.
4,000 K LED illumination (CRI: 92).
5X magnification.

Made in Germany.