

RIESTER ELITEVUE LED OTO-OPHTHALMO 3.5V - 2210 -203



Product Code	31863
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC	4045396209855
---------	---------------

Description

RIESTER ELITEVUE LED OTO-ophthalmo 3.5 V - 2210-203

Both with C-handle for rechargeable Li-Ion battery
EliteVue offers twice the field of view and 66% greater magnification (at 15 mm distance) compared to standard otoscopes.
Improved view of auditory canal and tympanic membrane guarantees more efficient diagnosis.
Suitable for throat illumination.
4,000 K LED illumination (CRI: 92).
5X magnification.

Supplied in a case with specula

Made in Germany.