

CRYO PNEUMATIC KNEE ORTHOSIS - universal



Product Code	42030
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2737096
CND	Y090618
EAN/UPC	8023279420302
PARAF	950080135

Description

CRYO PNEUMATIC KNEE ORT HOSIS

Universal cryo pneumatic knee orthosis indicated for following ACL repair, post total joint surgery, meniscal tears, patella realignment or during rehabilitation.

ROM hinges with malleable aluminium uprights offering full flexion and extension control with settings at 0°, 15°, 30°, 60° and 90°. Made of plush foam laminate with hook and loop closures, reusable gel packs feature soft fabric lining and nontoxic gel. Gel pack temperature: -20°C + 50°C.

Size fits up to 61 cm circumference.

Supplied with 2 ROM hinges with plastic covers and extension straps, 2 hot/cold gel packs, air compression body wrap, hand pump.